

E-SCOPE[®] Clinical Model



The **E-Scope (Electronic Stethoscope)**

Clinical Model allows medical professionals to auscultate with increased volume by amplifying sounds up to 30 times louder than an acoustic scope. The E-Scope is popular around the world in clinics, hospitals, research laboratories, and telemedicine settings. Additionally, the E-Scope can be paired with an accessory cord to broadcast sounds through speakers or other audio outputs.



One-handed control over volume adjustments



Specialist adult diaphragm



Automatic shut off after 1.5 to 2 minutes



Audio output jack



Battery operated (single AAA battery)



Supplied with accessory pack

Cardionics
The Heart of Auscultation

Phone: 281-488-5901, 800-364-5901 (USA)
Fax: 281-488-3195
Email: info@cardionics.com
Web: www.cardionics.com
Address: 910 Bay Star Blvd, Webster, TX 77598 USA

*E-Scope Catalog #718-7700 (Clinical model)
#718-7710 (Hearing impaired model)
Technical Support: Unlimited
Warranty: One year – parts and labor*