

E-SCOPE®

HEARING IMPAIRED



Used with hearing aid or hearing impairment



One-handed control over volume adjustments



Specialist adult diaphragm



Automatic shut off after 1.5 to 2 minutes



Audio output jack



Battery operated (single AAA battery)



Supplied with accessory pack



The **E-Scope Hearing Impaired Model** includes an assortment of auditory accessories, allowing users with hearing loss the ability to choose the appropriate listening device for their specific needs. This model is identical in functionality to the Clinical Model — amplifying sounds up to 30 times louder than an acoustic scope. Separate headphones and accessories are offered to ensure we find the right solution for any user with hearing loss.

E-Scope Catalog # 718-7700 (Clinical model)

718-7710 (Hearing impaired model)

Technical Support: Unlimited

Warranty: One year – parts and labor

Cardionics
™ The Heart of Auscultation

Phone: 281-488-5901, 800-364-5901 (USA)
Fax: 281-488-3195
Email: info@cardionics.com
Web: www.cardionics.com
Address: 910 Bay Star Blvd, Webster, TX 77598 USA