



The **PAT BASIC** trainer engages the technology of the SimScope[®] stethoscope to simulate pediatric heart, lung, and bowel sounds utilizing a large sounds library. Using the provided SimScope stethoscope, users are able to access a large variety of pathological conditions in order to customize scenarios. By placing the SimScope over the anatomically correct locations of PAT BASIC, the programmed sounds will be played through the scope, allowing for students to employ and enhance their diagnostic skills.

The optional SimScope WiFi upgrade enables wireless communication between the SimScope and PAT BASIC trainer, allowing users to change sounds remotely or to use in conjunction with a scenario-based simulation.

PAT BASIC Catalog # 718-8835 (Light skin) 718-8837 (Dark skin) Technical Support: Unlimited Warranty: One year – parts and labor



 Phone:
 281-488-5901, 800-364-5901 (USA)

 Fax:
 281-488-3195

 Email:
 info@cardionics.com

 Web:
 www.cardionics.com

 Address:
 910 Bay Star Blvd, Webster, TX 77598 USA

- Utilizes the Extensive Cardionics Proprietary Library of Heart, Lung, Bowel, and Bruit Sounds
- Listen at Anatomically Correct Auscultation Sites
- Customizable to Fit Specific Scenarios and Programs
- Portable and Lightweight
- Easy to Use
- Light and Dark Skin Models Available
- WiFi Upgrade Available Allowing for Remote Selection of Sound Scenarios