

# SAM BASIC<sup>®</sup>

## STUDENT AUSCULTATION MANIKIN



The **SAM BASIC** trainer engages the technology of the SimScope<sup>®</sup> stethoscope to simulate adult heart, lung, and bowel sounds utilizing a large sounds library. Using the provided SimScope stethoscope, users are able to access a large variety of pathological conditions in order to customize scenarios.

By placing the SimScope over the anatomically correct locations of SAM BASIC, the programmed sounds will be played through the scope, allowing for students to employ and enhance their diagnostic skills.

The optional SimScope WiFi upgrade enables wireless communication between the SimScope and SAM BASIC manikin, allowing users to change sounds remotely or to use in conjunction with a scenario-based simulation.

*SAM BASIC Catalog # 718-8830 (Light skin) 718-8832 (Dark skin)*  
*Technical Support: Unlimited*  
*Warranty: One year – parts and labor*

- Utilizes the extensive Cardionics proprietary library of heart, lung, bowel, and bruit sounds
- Listen at anatomically correct auscultation sites
- Customizable to fit specific scenarios and programs
- Portable and lightweight
- Easy to use
- Light and dark skin models available
- WiFi upgrade available allowing for remote selection of sound scenarios

**Cardionics**  
The Heart of Auscultation 

Phone: 281-488-5901, 800-364-5901 (USA)  
Fax: 281-488-3195  
Email: [info@cardionics.com](mailto:info@cardionics.com)  
Web: [www.cardionics.com](http://www.cardionics.com)  
Address: 910 Bay Star Blvd, Webster, TX 77598 USA